

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) For self-test

Product Pictures

➤ One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-test)

1 Test /Kit (English)

Specifications

Info. of the Test Kit and Export Packing Cartons

Product Name	Specifications	Size (cm)	Weight/Kit (g)	Kit Quantity/Carton	Size of Carton (cm)	Weight per Carton (kg)
One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	1 Test/Kit	7*1.8*13 cm	25.5 g	300 Kits/Carton	49.5*40.5*30 cm	8.6 kg/Carton
	5 Tests/Kit	7*5.3*13 cm	76.6 g	120 Kits/Carton	55.5*41*30 cm	10.44 kg/Carton

Export Packing Cartons

Packing carton for 5 T/kit
Size: 55.5*41*30 cm

Packing carton for 1 T/kit
Size: 49.5*40.5*30 cm

Certificates and Clinical Evaluation

Clinical Agreement Study

Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-Test)

Sensitivity	97.06%
Specificity	98.71%
Total Percent Agreement	98.13%

The clinical performance of One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) was evaluated by testing a total of 480 nasal swab samples. It is compared to the results of RT-PCR assays. Overall study results were shown in the tables below.

Total		BGI's RT-PCR kit		
		Positive	Negative	Subtotal
Getein's kit	Positive	165	4	169
	Negative	5	306	311
	Subtotal	170	310	480

Positive percent agreement (Diagnostic sensitivity) = $165 / (165 + 5) \times 100\% = 97.06\%$ (95% CI: 93.30%-98.74%)

Negative percent agreement (Diagnostic specificity) = $306 / (306 + 4) \times 100\% = 98.71\%$ (95% CI: 96.73%-99.50%)

Total percent agreement = $(165 + 306) / 480 \times 100\% = 98.13\%$ (95% CI: 96.48%-99.01%)

**CE Certificate of Getein One Step Test
for SARS-CoV-2 Antigen (Colloidal
Gold) (For Self-Test)**

CERTIFICATE

EC Certificate No. 1434-IVDD-447/2021

**EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices**

Polish Centre for Testing and Certification certifies
that manufactured by:

GETEIN Biotech, Inc.

Nanjing, ul. Bofu Road, Luhe District 9, China

in vitro diagnostic medical devices
for self-testing

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)

*Ref. codes: CG20615, CG206152, CG206153, CG206155, CG206156, CG206157, CG206158, CG206159,
CG2061510, CG2061512, CG2061515, CG2061520, CG2061525*

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 30.07.2021 to 27.05.2024

The date of issue of the Certificate: 30.07.2021

The date of the first issue of the Certificate: 30.07.2021

Issued under the Contract No. MD-66/2021
Application No: 142/2021
Certificate bears the qualified signature.
Warsaw, 30.07.2021
Module A1

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Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (For Self-Test) has been approved by Malaysia MDA

 **PIHAK BERHAKASA PIHAK PENERIMA**
Medical Device Authority
KEMENTERIAN KESIHATAN MALAYSIA
Ministry of Health Malaysia
Ara 8, Prime 8, Prime Avenue 5,
Blok 3/47, Persiaran Apec,
43000 Cyberjaya, Selangor
Malaysia.
Tel: (+603)8330 0300
Faks: (+603)8330 0303
Pusat Runding: www.mda.gov.my
Email: mda@mda.gov.my

Ref :
Date : 10 November 2021

Dear Sir/Madam,

CONDITIONAL APPROVAL FOR IMPORTATION AND DISTRIBUTION OF MEDICAL DEVICE (COVID-19 SELF TEST KIT)

With reference to the above, I wish to inform that the Authority grants your establishment a conditional approval for the importation and distribution of medical device as listed in **Appendix 1**.

2. Please be informed that the validity of this conditional approval is from **10/11/2021** to **10/11/2022** and is subject to the following:

- Your establishment shall ensure that the medical device under this conditional approval complies with safety and performance requirements as stipulated in Medical Device Act 2012 (Act 737);
- Your establishment shall adhere to the conditions as stipulated in **Appendix 2**;
- The use of COVID-19 self test kit shall be limited for screening purpose only and all test results need further confirmation using RT-PCR.

3. This conditional approval for importation and distribution of this medical device is an interim measure in response to the current public health need during COVID-19 pandemic. This letter shall not be used for the purpose of promoting or advertising of the product and it does not exempt you from abiding to any laws or requirements by any other authorities of Malaysia.

Thank you,

(AHMAD SHARIFF BIN HAMBALI)
Chief Executive
Medical Device Authority
Ministry of Health Malaysia

Ref :
Date : 10 November 2021

Medical Device Details

Appendix 1

Name of Medical Device	: One Step Test for SARS-Cov-2 Antigen (Colloidal Gold)
Brand/Model	: GP
Identifier	: Refer Attachment 1
Sample type	: Nasal swab
Intended Use	: One Step Test for SARS-Cov-2 Antigen (Colloidal Gold) is intended for the qualitative detection of SARS-Cov-2 antigens in human nasal swab samples. This test is used for individuals suspected of COVID-19 within the first seven days of the onset of symptoms, such as headache, fever, cough, sore throat, loss of the sense of smell or taste, shortness of breath, muscle pain. Meanwhile the test can also be applied for individuals without symptoms. Results are for the identification of SARS-Cov-2 antigen. Positive results indicate the presence of SARS-Cov-2 antigens, but individual history and other diagnostic information is necessary for determine infection status. Negative results do not rule out SARS-Cov-2 infection. Negative results for individuals with symptoms similar to COVID-19 infection for more than seven days should be treated as negative possibly. If necessary, it should be confirmed by molecular assay. One Step Test for SARS-Cov-2 Antigen (Colloidal Gold) intended to be used to help the diagnosis of SARS-Cov-2 infection in upper respiratory samples during the acute phase of infection
Brief Description	: One Step Test for SARS-Cov-2 Antigen (Colloidal Gold) is intended for the qualitative detection of SARS-Cov-2 antigens in human nasal swab samples. This test is used for individuals suspected of COVID-19 within the first seven days of the onset of symptoms, such as headache, fever, cough, sore throat, loss of the sense of smell or taste, shortness of breath, muscle pain. Meanwhile the test can also be applied for individuals without symptoms. Results are for the identification of SARS-Cov-2 antigen. Positive results indicate the presence of SARS-Cov-2 antigens, but individual history and other diagnostic information is necessary for determine infection status. Negative results do not rule out SARS-Cov-2 infection. Negative results for individuals with symptoms similar to COVID-19 infection for more than seven days should be treated as negative possibly. If necessary, it should be confirmed by molecular assay. One Step Test for SARS-Cov-2 Antigen (Colloidal Gold) intended to be used to help the diagnosis of SARS-Cov-2 infection. This test is used to help the diagnosis of SARS-Cov-2 infection.
Lot Number	: PSC213020W
Manufacturer's name	: Getein Biotech Inc. P.S.

**Getein One Step Test for SARS-CoV-2
Antigen (Colloidal Gold) (For Self-Test)
has been approved by Thailand FDA**

[illegible][illegible]

Getein One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) has passed the PEI Evaluation

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
**One Step Test for SARS-CoV-2 Antigen
(Colloidal Gold)**

Getein Biotech, Inc.

Getin One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) has passed the phase 3a validation of British DHSC

TopicsGovernment activity

Coronavirus (COVID-19) | Guidance and support

Home > Coronavirus (COVID-19) > Testing for coronavirus (COVID-19) > Assessment and procurement of coronavirus (COVID-19) tests

Department of Health & Social Care

Guidance

Outcome of the evaluation of rapid diagnostic assays for specific SARS-CoV-2 antigens (lateral flow devices)

Updated 11 October 2021

Background

Since its establishment in August 2020, the joint PHE Porton Down and University of Oxford SARS-CoV-2 lateral flow antigen test validation cell has evaluated over 140 lateral flow devices that have been referred by the Department of Health and Social Care (DHSC).

Approximately 30% of the tests that were referred for validation met the standards for phase 2 validation, which are set out in the [protocol](#). PHE Porton Down subsequently performed phase 3 testing to assess whether the lateral flow devices that passed phase 2 displayed performance characteristics desirable for mass population, community-based testing. The desirable performance characteristics are:

- very high specificity
- very high sensitivity against viral loads associated with infectiousness

The lateral flow devices that display the desirable performance characteristics are summarised in table 1 below.

Getin Biotech One Step Test for SARS-CoV-2 Antigen

Pass10 March 2021

CE Common List

EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY
Public health, country knowledge, crisis management
Health Security

EU health preparedness:

**A common list of COVID-19 rapid antigen tests,
including those of which their test results are mutually
recognised, and a common standardised set of data to
be included in COVID-19 test result certificates**

Agreed by the Health Security Committee

This document was agreed by the HSC on 17 February 2021

Annex I

Common list of rapid antigen tests

A first update was agreed by the HSC on 10 May 2021

A second update was agreed by the HSC on 16 June 2021

Annex II

Common standardised data set of to be included in COVID-19 test result certificates

An update to Annex II was agreed by the HSC on 19 March 2021

Getein Biotech, Inc.	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	Yes	97.06% sensitivity 98.71% specificity Nasal swab	DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 90% at <Ct30 and 100% at <Ct25)		DE ^[2]		DE ^[2]		Yes (2183)	No
----------------------	---	-----	--	--	--	-------------------	--	-------------------	--	---------------	----

Certificate of ISO 13485:2016

bsi.

By Royal Charter

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Getein Biotech, Inc.
No.9 Bofu Road
Luhe District
Nanjing
Jiangsu
211505
China

Facility ID Number: F004902

Holds Certificate No:

MDSAP 728434

Statement of Conformity: The company listed on this certificate has been audited to and found to conform with the following criteria: ISO 13485:2016 and Australia - Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure; Brasil - RDC ANVISA n. 16/2013, RDC ANVISA n. 23/2012, RDC ANVISA n. 67/2009; Canada - Medical Devices Regulations - Part 1 - SOR 98/282; Japan - MHLW Ministerial Ordinance 169, Article 4 to Article 68, PMD Act; USA - 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D

Please see scope page.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-10-18

Effective Date: 2020-10-18

Expiry Date: 2023-07-25

MDSAP

MEDICAL DEVICE SINGLE AUDIT PROGRAM
BSI Group America Inc. is an MDSAP authorized auditing organization

Page: 1 of 2

...making excellence a habit.™

Thanks for Your Attention